

User Instruction

December 17, 2025

Contents

1	Installation	2
2	Usage	3
3	Troubleshooting	3

1 Installation

To install the software, follow these steps:

APK Sideloading via SideQuest:

1. Download SideQuest from the official website <https://sidequestvr.com/>.
2. Install SideQuest on your computer by following the installation instructions provided on the website.
3. Connect the Meta Quest 3 headset to your computer using a USB-C cable.
4. Enable Developer Mode on your Meta Quest 3 headset by following the instructions on the Oculus website <https://developer.oculus.com/documentation/native/android/mobile-device-setup/>.
5. Open SideQuest and ensure that it recognizes your connected Meta Quest 3 headset.
6. Sideload the application `Quest_XueWei.apk` onto your Meta Quest 3 headset using SideQuest. You can find the APK file in the `.../project_source/` directory.
7. Once the installation is complete, disconnect the headset from your computer.

Server Setup:

1. Ensure you have Python 3.7 or higher installed on your server machine.
2. Make sure you are in the directory `.../project_source/` where the `server.py` and `server_requirements.txt` files are located.
3. Install the required Python packages by running the following command in your terminal:

```
pip install -r server_requirements.txt
```

4. Start the server by running the following command in your terminal:

```
nohup python3 -m uvicorn server:app \  
--host 0.0.0.0 \  
--port 8000 \  
> uvicorn.log 2>&1 &
```

5. Ensure that the server is running and accessible from your Meta Quest 3 headset.

Once the installation and server setup are complete, you can launch the application on your Meta Quest 3 headset and connect to the server to start using the software.

2 Usage

- To start using the application, simply upload your desired patient image to the server through the web interface provided by the application. The application will process the image and provide the necessary outputs based on the functionalities implemented.
- Upon completing the upload, log in to the Meta Quest 3 application. You should see the processed patient image displayed within the VR environment.
- To interact with the image, use the Meta Quest 3 controllers or hand tracking features to navigate, zoom, and manipulate the image as needed.
- Once a new image is uploaded to the server, first click “Return” to go back to the main menu, then click the refresh button to load the latest image.

3 Troubleshooting

If the initial loading screen persists for an extended period, please try the following steps:

1. Ensure that your Meta Quest 3 headset is properly connected to the server and that the server is running without issues.
2. Check your internet connection to ensure that there are no connectivity problems.
3. Restart the Meta Quest 3 application and try reconnecting to the server.
4. If the problem persists, consider restarting both the Meta Quest 3 headset and the server.

If hand tracking is not functioning correctly, please try the following steps:

1. Ensure that hand tracking is enabled in the Meta Quest 3 settings.
2. Make sure that your hands are within the field of view of the headset’s cameras.
3. Make sure that you are in a well-lit environment to improve hand tracking accuracy.
4. Make sure that your background is not too cluttered, as this can interfere with hand tracking.
5. Restart the Meta Quest 3 application to refresh the hand tracking functionality.